




SYM Knife®

SYM Quadra / Tres

for Microtome & Ultra Microtome



**Huge boat collects
numerous sections at once**



Specialized specifications for collecting serial sections

Features of SYM Quadra

Point 1: 3-inch wafers can be installed as is without cutting




The structure inside the boat has the suitable form for setting wafers, and a 3-inch silicon wafer fits perfectly.

Point 2: A lot of ribbons can be placed in one piece



As shown in the photos above, many serial sections can be gathered on a wafer.

Point 3: It has a drainage mechanism for collection



Serial sections can be collected easily, wrinkle-free, and quickly using a drainage mechanism.

The Line-Up of SYM Quadra

Item No.	Width	Angle
QT2035	2.0mm	35°
QT2045	2.0mm	45°
QT3035	3.0mm	35°
QT3045	3.0mm	45°
QT4035	4.0mm	35°
QT4045	4.0mm	45°

Make time your friend!



SYNTEK CO., LTD (Japan)

