






# SYM Knife Multi Grid Holder

## SYM Knife

### This is a diamond knife what I wanted !

(Accessory function holder for collecting samples)



There're some more advantages.

Including some of Universal design

- Everybody uses it in the same way
- No matter who you're right or left-handed
- Intuitively know how to use it

1

Normal method	In Multi Grid Holder
Likely to get wrinkled	Hardly becomes wrinkled

Crushed!

2

Conventional way	Multi Grid Holder
Stick out-through support membrane	Swoop-swoop-slide

Thump!

Unlikely to break support membrane

Pushing sample on the grid

3

You can collect just needed ones from many sections

It's easy to put a sample at the target position because grid is fixed. Also, you collect them on the same direction.

4

We uploaded our MGH video. Check it out for details.

Click "Multi grid holder" on YouTube

5



SYNTEK CO., LTD.

info@syntek.co.jp  
<http://www.syntek.co.jp/en>

Patent-Pending